

# NOBILTRAY

— Description

NOBILTRAY is a preformed light curing tray indicated for making individual tray and can also work as temporary base plate in dental labs. As light curing material, it can be fast cured in 3 minutes in a unit. The product is characterized by easy operation, high precision, low deformation rate, high plasticity and pleasant mint flavour.



— Application

- Making individual tray



Types Features	Light Curing Tray	Self-cured resin Tray	Metal Tray
Operation	preformed → ready for use and time-saving	need mixing and kneading first → complicated and time-consuming	-
Plasticity	soft but with certain strength → easy to trim and shape when opening windows	too soft to work with in dough stage → difficult to trim and shape when opening windows	hard and rigid metal material → difficult to open window
Accuracy	high hardness after curing → low deformation rate ensuring accuracy	existence of contraction → unstable and inaccurate	specific and uniform size → fail to match different oral forms and inaccurate
Smell	mint flavour → pleasant smell	residual stimulative monomer → pungent smell	-

- Making temporary base plate

Types Features	Light Curing Base Plate	Wax base Plate
Operation	preformed → ready for use and time-saving	need heating first → time-consuming
Accuracy	high hardness after curing → low deformation rate ensuring accuracy	vulnerable to temperature → unstable and inaccurate



— Main Advantages

- Preformed → easy handling
- 3-minute fast curing → quick curing
- High hardness and fit different oral forms → accurate
- Soft → high plasticity
- Mint flavour → pleasant smell
- Various choices → pink or blue, thicker or thinner, all available

— Specification

Standard package



50pcs/box

Sample package



2pcs/bag



Follow @VinciSmile

VINCISMILE GROUP LLC

Add: 20524 Carrey Rd. Walnut, CA 91789, USA  
Email: [marketing@vincismile.com](mailto:marketing@vincismile.com)  
Web: [www.vincismile.com](http://www.vincismile.com)