

i-Vinci SCANNER

Rizhao | Beijing | Shanghai | Qingdao | LA

Manufactured for:

VINCISMILE GROUP LLC

Add: 20524 Carrey Rd. Walnut, CA 91789, USA

Email: marketing@vincismile.com

Web: www.vincismile.com

Manufacturer:

Ziyang Freqty Medical Equipment Co., Ltd.

Add: Zone A, B, C, 4 Floor, Building A, No.3, Xiandai Road,
Ziyang City, Sichuan Province, China

i-Vinci Smart

Open A New Era of Digital Dentistry

Only **138g** 221x27x25mm



SPECIFICATION

Technology	Staring Scanning with Coded Structure Light
Anti Fog System	Intelligent Heating
Accuracy	7 μ m
Precision	5 μ m
Handpiece Weight	138g (192g with cable)
Handpiece Dimension	221x27x25mm
Handpiece Enclosure	Aviation Aluminum Alloy
Scan Field	18x16mm
Scan Depth	0~20mm
True Color	Full HD
Tips	x6 (Autoclave)
Tip Autoclave Cycle	50~100 times
Scanning Control	Button
Motion Sensing	Support
Data Interface	USB3.0
Cable Length (m)	1.8m

Smaller Tip

Faster Scan

Higher Accuracy

Motion Sensing

Super Quiet

